



Recombinant Human Apolipoprotein O-like (APOOL)

Product Code	CSB-BP001949HU
Abbreviation	APOOL
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6UXV4
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	AKQEE SKKQLVKPEQ LPIYTAPPLQ SKYVEEQPGH LQMGFASIRT ATGCYIGWCK GYVVFVKNGI MDTVQFGKDA YVYLKNPPRD FLPKMGVITV SGLAGLVSAR KGSKFKKITY PLGLATLGAT VCYPVQSVII AKVTAKKVYA TSQQIFGAVK SLWTKSSKEE SLPKPKEKTK LGSSSEIEVP AKTTHVLKHS VPLPTELSSE AKTKSESTSG ATQFMPDPKL MDHGQSHPED IDMYSTRS
Source	Baculovirus
Target Names	APOOL
Protein Names	Recommended name: Apolipoprotein O-like Alternative name(s): Protein FAM121A
Expression Region	26-268
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.