



Recombinant Human Zinc finger protein 562 (ZNF562)

Product Code	CSB-BP026843HU
Abbreviation	ZNF562
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6V9R5
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MSAFDMSHGF FPREPICPFE EGTKIGTMVE DHRSNSYQDS VTFDDVAVEF TPEEWALLDT TQKYLRYDVM LENYMNLASV DFFFCLTSEW EIQPRTKRSS LQQGFLKNQI FTGIQMQRTR YSGWKLCENC GEVFSEQFCL KTHMRAQNGG NTFEGNCYGK DSISVHKEAS IGQELSKFNP CGKVFTLTPG LAVHLEILNG RQPYKCKEKG KGFKYFASLD NHMGIHIGEK LCEFQECERA ITTSSHLKQC VAVHTGKKSE KTKNCGKSFT NFSQLSAHAK THKGEKSFEK KECGRSFRNS SSFNVHIQIH TGIKPHKCTE CGKAFTRSTH LTQHVTRHTG IKPYECKEKG QAFTQYTGLA IHIRNHTGEK PYQCKECKGA FNRSSTLTQH RRIHTGEKPY ECVECGKTFI TSSHRSKHLK THSGER
Source	Baculovirus
Target Names	ZNF562
Protein Names	Recommended name: Zinc finger protein 562
Expression Region	1-426
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.