



Recombinant Human Zinc finger protein 367 (ZNF367)

Product Code	CSB-EP026703HU-B
Abbreviation	ZNF367
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7RTV3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MIRGFEAPMA ENPPPPPPV IFCHDSPKRV LVSVIRTTPI KPTCGGGGEP EPPPPLIPTS PGFSDFMVYP WRWGENAHNV TLSPGAAGAA ASAALPAAAA AEHSGLRGRG APPPAASASA AASGGEDEEE ASSPDSGHLK DGIRRGPRRA DTVRDLIN EG EHSSSRIRC N ICNRVFPREK SLQAHKRTH GERPYLCDYP DCGKAFVQSG QLKTHQRLHT GEKPFVCS EN GCLSRFTHAN RHCPKHPYAR LKREEPTDTL SKHQAADNKA AAEWLARYWE MREQRTPTLK GKL VQKADQE QQDPLEYLQS DEEDDEKRG A QRRLQEQRER LHGALALIEL ANLTGAPLRQ
Source	E.coli
Target Names	ZNF367
Protein Names	Recommended name: Zinc finger protein 367 Alternative name(s): C2H2 zinc finger protein ZFF29
Expression Region	1-350
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.