



Recombinant Mouse Inositol polyphosphate multikinase (Ipmk)

Product Code	CSB-BP011774MO
Abbreviation	Ipmk
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7TT16
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	AAEPPALRL RPPGSTGDSP PVPRLLGGCV PLSHQVAGHM YGKDKVIGILQ HPDGTVLKQL QPPPRGPREL EFYTMVYAAD CADAVLLELR KHLPKYYGVW SPPTAPNDVY LKLEDVTHKF NKPCIMDVKI GRKSYDPFAS SEKIQQQVSK YPLMEEIGFL VLGMRVYHLH SDSYETQNH YGRGLTKETL KEGVSKFFHN GFCLRKDAIA ASIQKVEKIL QWFENQKQLN FYASSLLFVY EGSSQPATTK ANDRTLGRF LSKGPLTDAD GLECNFFFHL FGAPPNGMSV GKSLSKAYSR HRKLYAKKHQ SQTSLKVETL EQDNGWRSMS QEHLNGNVLA QLEKVFYHLP AGRPEIPEAE VRMIDFAHVF PSNTVDEGYV YGLKHIAVL RSILDS
Source	Baculovirus
Target Names	Ipmk
Protein Names	Recommended name: Inositol polyphosphate multikinase EC= 2.7.1.151 Alternative name(s): Inositol 1,3,4,6-tetrakisphosphate 5-kinase
Expression Region	2-396
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.