



# Recombinant Human Zinc finger protein 544 (ZNF544)

<b>Product Code</b>	CSB-EP026826HU
<b>Abbreviation</b>	ZNF544
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q6NX49
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MEARSMLVPP QASVCFEDVA MAFTQEEWEQ LDLAQRTLIR EVTLETWEHI VSLGLFLSKS DVISQLEQEE DLCRAEQEAP RDWKATLEEN RLNSEKDRAR EELSHHVEVY RSGPEEPPSL VLGKVQDQSN QLREHQENSL RFMVLTSERL FAQREHCELE LGGGYSLPST LLLPTTLPT STGFPPKNSQ VKELKQNSAF INHEKNGADG KHCESHQCAR AFCQSIYLSK LGNVETGKKN PYEYIVSGDS LNYGSSLCFH GRTFSVKKSD DCKDYGNLFS HSVSLNEQKP VHFQKSYEC DECRETCSES LCLVQTERSG PGETPFRCEE RCAAFFPMAS FSDCNIIQTT EKPSVCNQCQ KSFSCCKLIH QRTHTGEKPF ECTQCGKSFS QSYDLVIHQR THTGEKPYEC DLGKSFTQR SKLITHQRIH TGEKPYQCIE CRKSFRWNSN LIVHQRHTG EKPYECHCG KSFSQSYELV THKRHTGEEK PFKCTQCGKS FSQKYDLVVH QRTHTGEKPY ECNLCGKSFS QSSKLITHQR IHTGEKPYQC IECGKSFRWN SNLVIHQRIH TGEKPYDCTH CGKSFSQSYQ LVAHKRHTGT EKPYECEG KAFNRSTQLI RHLQIHTGEEK PYKCNQCNKA FARSSYLVMH QRTHTGEEKPF ECSQCGKAFS GSSNLLSHHR IHSGEKPYEC SDCGKSFRQQ SQLVVHRRTH TGEKP
<b>Source</b>	E.coli
<b>Target Names</b>	ZNF544
<b>Protein Names</b>	Recommended name: Zinc finger protein 544
<b>Expression Region</b>	1-715
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.