



Recombinant Danio rerio Cerebellar degeneration-related protein 2-like (cdr2l)

Product Code	CSB-MP005104DIL
Abbreviation	cdr2l
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6NZT2
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MLRAGR MDEF VTEED EPWYD QRDLEQDLHL AAELGKTLLE RNKELEDSLQ QMYINNEEQV QEIEYLTKQM EMLREMNEQH AKVYEQLDVT ARELEITNEK LVLESKGSQQ KIDRLTGTME MLQGQVDTLT ARVEELRTLE ELRVRREKKE RRKTVHSFPC LKELCTAPRY EDSFMVSDPG SGDLEECQPA DEENEHLRVM VSSLRAAVAA ERGRREAAER ECAAVLQEFE RLEQRLLGAE SCQRRVHELE AELQEMQQLR KSRVCLLSGE EGGLEQTL LNC APETDTPDEA VMAQNGDAND GAAPVRKSCS DTALDAISAV DASGRRKGSY ALHANGVRKR GMSILREVDE QYHALLEKYE ELLGKCRHE ESLRHAEVQT SRPVS RDPSM KECRVAEPQQ PPPTPPQTPS TPEALEGISR QVEAVDKRLS QNTPEYKALF KEIFSRLQRT KSDINSTKGR KGGK
Source	Mammalian cell
Target Names	cdr2l
Protein Names	Recommended name: Cerebellar degeneration-related protein 2-like
Expression Region	1-464
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.