



Recombinant Human Zinc finger protein 549 (ZNF549)

Product Code	CSB-BP026830HU
Abbreviation	ZNF549
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6P9A3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAEAALVITP QIPMVTEEFV KPSQGHVTFE DIAVYFSQEE WGLLDEAQRCLYHDVMLENF SLMASVGCLH GIEAEEAPSE QTLAQGVVSQ ARTPKLGPSIPNAHSCEMCI LVMKDILYLS EHQTLPWQK PYTSVASGKW FSFGSNLQQHQNQDSGEKHI RKEESSALLL NSCKIPLSDN LFPCKDVEKD FPTILGLLQHQTTHSRQEYA HRSRETFQQR RYKCEQVFNE KVHVTEHQRVHTGEKAYKRR EYGKSLNSKY LFVEHQRTHN AEKPYVCNIC GKSFLHKQTLVGHQQRIHTR ERSYVCIECG KSLSSKYSLV EHQRTHNGEK PYVCNVCGKSFRHKQTFVGH QQRIHTGERP YVCEMGKSF IHSYDRIRHQ RVHTGEGAYQCSECGKFSFY KQSLLDHHRI HTGERPYECK ECGKAFIHKK RLEHQRIHTGEKPYVCIIC GKSFIRSSDY MRHQRIHTGE RAYECSDCGK AFISKQTLHKHHKIHTRERP YECSECGKGF YLEVKLLQHQ RIHTREQLCE CNECGKVFESHQKRLLEHQKV HTGEKPCECS ECGKCFRHRT SLIQHQKVHS GERPYNCTACEKAFIYKNKL VEHQRIHTGE KPYECGKCGK AFNKRYSLVR HQKVHTEEP
Source	Baculovirus
Target Names	ZNF549
Protein Names	Recommended name: Zinc finger protein 549
Expression Region	1-640
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.