



# Recombinant Human Suprabasin (SBSN)

<b>Product Code</b>	CSB-BP020734HU
<b>Abbreviation</b>	SBSN
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q6UWP8
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	DDPIE KVIEGINRGL SNAEREVGKA LDGINSGITH AGREVEKVFN GLSNMGSHGTG KELDKGVQGL NHGMDKVAHE INHGIGQAGK EAEKLGHGTVN NAAGQVGKEA DKLIHHGVHHA GANQAGSEAG KFGQGVNDAA GQAGNEAGRF GQGVHHAAGQ AGNEAGRFGQ GVHHAAGQAG NEAGRFGQGA HHGLSEGWKE TEKFGQGIHH AAGQVGKEAE KFGQGAHHA GQAGNEAGRF GQGVHHGLSE GWKETEKFGQ GVHHTAGQVG KEAEKFGQGA HHAAGQAGNE AGRFGQGAHH AAGQAGNEAG RFGQGVHHGL SEGWKETEFK GQGVHHAASQ FGKETEKLGH GVHHGVNEAW KEAEKFGQGV HHAASQVGKE EDRVVQGLHH GVSQAGREAG QFGHDIHHTA GQAGKEGDIA VHGVPQGVHE AGKEAGQFGQ GVHHTLEQAG KEADKAVQGF HTGVHQAGKE AEKLGQGVNH AADQAGKEVE KLGQGAHHA GQAGKELQNA HNGVNVQASKE ANQLLNGNHQ SGSSSHQGGG TTTPLASGAS VNTPFINLPA LWRSVANIMP
<b>Source</b>	Baculovirus
<b>Target Names</b>	SBSN
<b>Protein Names</b>	Recommended name: Suprabasin
<b>Expression Region</b>	26-590
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.