



Recombinant *Zalophus californianus* Apolipoprotein E (APOE)

Product Code	CSB-EP001936ZAC
Abbreviation	APOE
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7M2U7
Product Type	Recombinant Protein
Immunogen Species	<i>Zalophus californianus</i> (California sealion)
Purity	>85% (SDS-PAGE)
Sequence	DV EPESPLEENL EPELEPKREL EQEVEPEAGW QAGQPWELAL ARFWDYLRWV QTLSDQVQEE VLSNQVTQEL TTLMEETMKE IKAYRAELEE QLGPMASETQ ARVAKELQAA QARLRSDMED VRTRLSQYRG EVQAMLGQST EELRARFASH MRKLRKRVLR DAEDLQKRLA VYRAGVREGA ERSVSTIRER LWPLLEQART RHAKVDALAT QPLRERVNAL GQQLRGRLEE VGSRRARSHLD EVREQMEEVQ AKMEEQANQM RQQAFAFQAR LKGWFEPLVE DMQRQWAVLV EKVQAAVGTS PTTPPVETK
Source	E.coli
Target Names	APOE
Protein Names	Recommended name: Apolipoprotein E Short name= Apo-E
Expression Region	19-329
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.