



Recombinant Chromobacterium violaceum Uracil-DNA glycosylase (ung)

Product Code	CSB-EP025641CKA-B
Abbreviation	ung
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7P111
Product Type	Recombinant Protein
Immunogen Species	Chromobacterium violaceum (strain ATCC 12472 / DSM 30191 / JCM 1249 / NBRC 12614 / NCIMB 9131 / NCTC 9757)
Purity	>85% (SDS-PAGE)
Sequence	MTYLAPALAR IHPDWEAVLD QAAVAARLAE IDRLLVLQQQ DGKTLFPPPP QVFNALAFAA PADVKVVILG QDPYHGDGEA MGLSFSVPDG ARVPPSLRNI YKELAADLGL GVPASGDLTH WAQQGVLLLN SVLTVERDKA GSHGKLGWQT VSDALIDAVN RDNPGCVFLL WGNWAQTKAE RIDASRHLVL TAAHPSPLSA SRGFHGCRRHF SQVNAWLIAR GRQPVRWATA PAAQSCLF
Source	E.coli
Target Names	ung
Protein Names	Recommended name: Uracil-DNA glycosylase Short name= UDG EC= 3.2.2.27
Expression Region	1-238
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.