



Recombinant Danio rerio PHD finger protein 10 (phf10)

Product Code	CSB-EP017898DIL-B
Abbreviation	phf10
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6NWE1
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MAAVLSRPC DSNPATPGTQ SIKDDIEENS SDGSQAPKRR RMGSGDSSRS CDTSTHELGP TYFPAENLTE YKWPPDDTGE YYMLQEQVSE YLGVTSFKRK YPEMERRDLS HKEKLYLREQ NVITETQCTL GLTALRSDEV IDLMIKEYPA KHAEYSVILQ ERERQRITKE YNVSTLSVQM QQQNPQKVEA SKVPEYIKKA AKKAAEFNSN FNRERMEERR AYFDLQTHII QVPQGKYKIL PPERTKTGPF PVALIPGQFQ EYKRYSPNE LRYLPLNTAL FEPPLDPELP ALSDGSDSD GDDGKEDGKK NKGSSDSSSG NTSDGDSQDS GIHSKVKSKD RPGAQKDGDT PRSIQHKSVP GYKPKVIPNA ICGICQKGKE ANKRGKPEAL IHCSQCQNSG HPSCLDMSVD LVAKIKMYPW QCMECKTCTV CQQPHHEEEM MFCDKCDRGF HTFCVGMDSI PMGCWVCDLC SKDISTPQKK GQTKTPKKAK
Source	E.coli
Target Names	phf10
Protein Names	Recommended name: PHD finger protein 10
Expression Region	1-490
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.