



# Recombinant Danio rerio Phospholipase A1 member A (pla1a)

<b>Product Code</b>	CSB-YP018084DIL
<b>Abbreviation</b>	pla1a
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q6NYZ4
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Danio rerio (Zebrafish) (Brachydanio rerio)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	DDEEASASK CADFNNTTWL EYRQATKLQV QYLLLTRKNA NCASLFTQDC LNHTQKHTAY FNSSLPTKVI VHGYRALGSK PSWVSGLAQA LLREEDVNVL VVDWVYGASF AYNLVVENYK EVAVQISVLI NQLTKYGSTL ESFHFIGVSL GAHVSGFVGT LFHGKLGKIT GLDPAGPMFK SADPFDRDLS SDALFVEAIH TDSDFYFGISI PVGHVDFFLN GGMDQAGCAR SRFASMYGYV ICDHMRALHV YMSALNGSCP LIGFPCSGYE EFLAGKCITC DDPFNGTCPQ IGLLKNSGIT ATPLPNQEKV YLLTTASGPF CAHHILVELN VSRLDKTAEV QLILKSIGHP ETELNLKLYT DETRYKTVAH HPEQLCKIDS MQLINTGGRF YRQGDHFEY ICISEIPQTR GMDPLCVKNI HIGRGPVWSH DFVQVC
<b>Source</b>	Yeast
<b>Target Names</b>	pla1a
<b>Protein Names</b>	Recommended name: Phospholipase A1 member A EC= 3.1.1.-
<b>Expression Region</b>	22-456
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.