



Recombinant Human Arylacetamide deacetylase-like 2 (AADACL2)

Product Code	CSB-BP001016HU
Abbreviation	AADACL2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6P093
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	FY TPMPDNIEES WKIMALDAIA KTCTFTAMCF ENMRIMRYEE FISMIFRLDY TQPLSDEYIT VTDTFVDIP VRLYLPKRKS ETRRRAVIYF HGGGFCFGSS KQRAFDFLNR WTANTLDAVV VGVDYRLAPQ HHFPAQFEDG LAAVKFFLLE KILTKYGVDP TRICIAGDSS GGNLATAVTQ QVQNDAEIKH KIKMQVLLYP GLQITDSYLP SHRENEHGIV LTRDVAIKLV SLYFTKDEAL PWAMRRNQHM PLESRHLFKF VNWSILLPEK YRKDYVYTEP ILGGLSYSPL GLTDSRALPL LANDSQLQNL PLTYILTCQH DLLRDDGLMY VTRLRNVGVQ VVHEHIEDGI HGALSFMTSP FYLRLGLRIR DMYVSWLDKN L
Source	Baculovirus
Target Names	AADACL2
Protein Names	Recommended name: Arylacetamide deacetylase-like 2 EC= 3.1.1.-
Expression Region	19-401
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.