



Recombinant Human Cysteine-rich protein 3 (CRIP3)

Product Code	CSB-EP005971HU-B
Abbreviation	CRIP3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6Q6R5
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MSWTCPRCQQ PVFFAEKVSS LGKNWHRFCL KCERCHSILS PGGHAEHNGR PYCHKPCYGA LFGPRGVNIG GVGSYLYNPP TPSPGCTTPL SPSSFSPPRP RTGLPQGKKS PPHMKFTTGE TSLCPGCGEP VYFAEKVMSL GRNWHRPCLR CQRCHKTLLTA GSHAEHDGVP YCHVPCYGYL FGPKGGQPHP RHWDGMYMPE VWHVHGLWVC VDNFPCG
Source	E.coli
Target Names	CRIP3
Protein Names	Recommended name: Cysteine-rich protein 3 Short name= CRP-3 Alternative name(s): Chromosome 6 LIM domain only protein Short name= h6LIMo
Expression Region	1-217
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.