



Recombinant Mouse Forkhead box protein O6 (Foxo6)

Product Code	CSB-EP008840MO
Abbreviation	Foxo6
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q70KY4
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	<p>MAAKLRAHQV DVDPDFAPQS RPRSCTWPLP QPDLAGDEDG ALGAGVAEGS EDCGPERRAT APAMAPAPPL GAEVGPLRKA KSSRRNAWGN LSYADLITKA IESAPDKRLT LSQIYDWMVR YVPYFKDKGD SNSSAGWKNS IRHNLSLHTR FIRVQNEGTG KSSWWMLNPE GGKTGKTPRR RAVSMDNGAK FLRIK GKASK KKQLHLPERS PDDSPPGAPV PGPLSASAKW AASPASHASD DYEAWADFRG SRRPLLGEAA ELEDDEALEA LAPSSPLMYP SPASALSPAL GARCPGELPR LAELGGPLGL HGGGVAGLPD ALLDGAQDAY GPRARAGTPS YFGSCKASAY GGGGGFGPPA LGSLRRLPMQ TIQENKQASF VQAAAPFRPG ALPALLPPPP PAPERPGLLG APGELALAGA AAAYPGKGAA PYAPPAPSR S ALAHPISLMT LPGEAGAAGL APPAHAAAFG GPPGGLLLLDA LPGPYAAAAA GPLGAGPDRF PADLDLDMFS GSLECDVESI ILNDFMDSDE MDFNFDSALP PPPPGLAGAP PPNQSWVPG</p>
Source	E.coli
Target Names	Foxo6
Protein Names	Recommended name: Forkhead box protein O6
Expression Region	1-559
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.