



# Recombinant Mouse Forkhead box protein O6 (Foxo6)

<b>Product Code</b>	CSB-MP008840MO
<b>Abbreviation</b>	Foxo6
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q70KY4
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MAAKLRAHQV DVDPDFAPQS RPRSCTWPLP QPDLAGDEDG ALGAGVAEGS EDCGPERRAT APAMAPAPPL GAEVGPLRKA KSSRRNAWGN LSYADLITKA IESAPDKRLT LSQIYDWMVR YVPYFKDKGD SNSSAGWKNS IRHNLSLHTR FIRVQNEGTG KSSWWMLNPE GGKTGKTPRR RAVSMDNGAK FLRIK GKASK KKQLHLPERS PDDSPPGAPV PGPLSASAKW AASPASHASD DYEAWADFRG SRRPLLGEAA ELEDDEALEA LAPSSPLMYP SPASALSPAL GARCPGELPR LAELGGPLGL HGGGVAGLPD ALLDGAQDAY GPRARAGTPS YFGSCKASAY GGGGGFGPPA LGSLRRLPMQ TIQENKQASF VQAAAPFRPG ALPALLPPPP PAPERPGLLG APGELALAGA AAAYPGKGAA PYAPPAPSR S ALAHPISLMT LPGEAGAAGL APPAHAAAFG GPPGGLLLLDA LPGPYAAAAA GPLGAGPDRF PADLDLDMFS GSLECDVESI ILNDFMDSDE MDFNFDSALP PPPPGLAGAP PPNQSWVPG
<b>Source</b>	Mammalian cell
<b>Target Names</b>	Foxo6
<b>Protein Names</b>	Recommended name: Forkhead box protein O6
<b>Expression Region</b>	1-559
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.