



Recombinant *Drosophila melanogaster* WAS protein family homolog 1 (wash)

Product Code	CSB-MP025972DLU
Abbreviation	wash
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7JW27
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MEESPYLHSP YQVAIIATDL HHEDTIIQAA QSLDCLHKTI NSIFERIDAR LARNGSKVED INNRVKRAQA KIDALVGSKR AIQIFAPARF PASDVLAPLP ATFPQVAANP LMEQQVDQLP QGTYSSSHSA DQKPDDADIF FHVRGDREQE SPLVAERKIT NRTAGLGILP AGGVRVPSL MRFNTNEFAY GEDLNAWKRS LPPQNARRVA SQSTQLTGEK QLAPAPHSLA HGTTKLATPA GDLRYNPAAL AAPAIDVPLD LPDLPGIAND LQYEPVEEQT PIAPSQQFGD LPELPDLGLE EQDIIVQAIA AQTHIPGPVR RKSVGQCPSV VTAAPPPPPP PPPPPPPPPP AQ TSAIPSP PFPTKGAVKP LPSLATPLN MPQPPATED PRSELMAAIR NAGGVHGGRL RSPAAAPLDV VDNSRSKAGG AVTGDLMADL HNKLMLRRKG ISGSQNPVEA TAGNPLMQQL SRVIPPPVQP RKGSKSSDEH SEDDGWN
Source	Mammalian cell
Target Names	wash
Protein Names	Recommended name: WAS protein family homolog 1 Alternative name(s): Protein washout
Expression Region	1-499
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.