



Recombinant Human Zinc finger protein 600 (ZNF600)

Product Code	CSB-YP026880HU
Abbreviation	ZNF600
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6ZNG1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MMKEVLSTGQ GNTEVIHTGT LQRYQSYHIG DFCFQEIEKE IHDIEFQCQE DERNGHEAPM TKIKKLTGST DQHDHRHAGN KPIKDQLGSS FYSHLPELHI IQIKGKIGNQ FEKSTSDAPS VSTSQRISPR PQIHISNNYG NNSPNSLLP QKQEVYMRK SFQCNEGKA FNCSSLLRKH QIPHLGDKQY KCDVCGKLFN HKQYLTCR CHTGEKPYKC NECGKSFSQV SSLTCHRRHLH TAVKSHKCNE CGKIFGQNSA LVIHKAHTG EKPYKCNECD KAFNQQSNLA RHRRIHTGEK PYKCECDKV FSRKSTLESH KRIHTGEKPY KCKVCDTAFT WNSQLARHKR IHTGEKTYKC NECGKTFCHK SSLVCHHRLH GGEKSYKCKV CDKAFWNSH LVRHTRIHSG GKPYKCNECG KTFGQNSDLL IHKSIHTGEQ PYKYEECEKV FSCGSTLETH KIIHTGEKPY KCKVCDKAFA CHSYLAKHTR IHSGEKPYKC NECKTFRLR SYLASHRRVH SGEKPYKCNE CSKTFQRSY LHCHRRHLHSG EKPYKCNECG KTFSHKPSLV HHRRLHTGEK SYKCTVCDKA FVRNSYLARH TRIHTAEKPY KCNECGKAFN QQSLSLHHR IHAGEKLYKC ETCDKVFSRK SHLKRHRRIH PGKKPYKCKV CDKTFGSDSH LKQHTGLHTG EKPYKCNECG KAFSKQSTLI HHQAVHGVGK LD
Source	Yeast
Target Names	ZNF600
Protein Names	Recommended name: Zinc finger protein 600
Expression Region	1-722
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.