



Recombinant Rat Heparanase (Hpse)

Product Code	CSB-BP010716RA
Abbreviation	Hpse
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q71RP1
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥85% (SDS-PAGE)
Sequence	KD VVDLEFYTKR LFQSVSPSFL SITIDASLAT DPRFLTFLGS PRLRALARGL SPAYLRFGGT KTDFLIFDPN KE
Source	Baculovirus
Target Names	Hpse
Protein Names	Recommended name: Heparanase EC= 3.2.-.- Alternative name(s): Endo-glucuronidase Cleaved into the following 2 chains: 1. Heparanase 8 kDa subunit 2. Heparanase 50 kDa subunit
Expression Region	29-102
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.