



Recombinant Human Fc receptor-like A (FCRLA)

Product Code	CSB-BP008560HU
Abbreviation	FCRLA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7L513
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	ASF ETLQCEGPVC TEESSCHTED DLTDAREAGF QVKAYTFSEP FHLIVSYDWL ILQGPAKPVF EGDLLVLRQC AWQDWPLTQV TFYRDGSALG PPGNREFSI TVVQKADSGH YHCSGIFQSP GPGIPETASV VAITVQELFP APILRAVPSA EPQAGSPMTL SCQTKLPLQR SAARLLFSFY KDGRIVQSRG LSSEFQIPTA SEDHSGSYWC EAATEDNQVW KQSPQLEIRV QGASSAAPP TLNPAPQKSA APGTAPPEAP GPLPPPPTPS SEDPGFSSPL GMPDPHLYHQ MGLLLKHMQD VRVLLGHLLM ELRELSGHRK PGTTKATAE
Source	Baculovirus
Target Names	FCRLA
Protein Names	Recommended name: Fc receptor-like A Alternative name(s): Fc receptor homolog expressed in B-cells Fc receptor-like and mucin-like protein 1 Fc receptor-like protein Fc receptor-related protein X Short name= FcRX
Expression Region	28-359
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.