



Recombinant Human Zinc finger protein 621 (ZNF621)

Product Code	CSB-BP026899HU
Abbreviation	ZNF621
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6ZSS3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MLQTTWPQES VTFEDVAVYF TQNQWASLDP AQRALYGEVM LENYANVASL VAFPPFKPAL ISHLERGEAP WGPDPWDTEI LRGISQGGES WIKNEGLVIK QEASEETELH RMPVGGLLRN VSQHFDFKRK ALKQTFNLNP NLILRGGMKF YECKEKGKIF RYNSKLIRHQ MSHTGEKPKF CKECGKAFKS SYDCIVHEKN HIGEGPYECK ECGKGLSNT ALTQHQRHT GEKPYECKEC GKAFRRSAAY LQHQLRHTGE KLYKCKECHKW AFGCRSLFIV HQRHTGEKP YQCKEKGKAF TQKIASIQHQ RVHTGEKPYE CKVCGKAFKW YGSFVQHKL HPVEKKPKV LGPSLVSPQC SSPAIPVLL QGSCSASAVA VPSLTFPHAV LIPTSGNFFM LLPTSGIPSS SAQIVRVFQG LTPTVKPSPV ILTPSSHSS
Source	Baculovirus
Target Names	ZNF621
Protein Names	Recommended name: Zinc finger protein 621
Expression Region	1-439
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.