



Recombinant Human RING finger protein 44 (RNF44)

Product Code	CSB-MP019893HU
Abbreviation	RNF44
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7L0R7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MRPWALAVTR WPPSAPVGQR RFSAGPGSTP GQLWGSPGLE GPLASPPARD ERLPSQQPPS RPPHLPVEER RASAPAGGSP RMLHPATQQS PFMVDLHEQV HQGPVPLSYT VTTVTTQGFP LPTGQHIPGC SAQQLPACSV MFSGQHYPLC CLPPPLIQAC TMQQLPVPYQ AYPHLISSDH YILHPPPPAP PPQPTHMAPL GQFVSLQTQH PRMPLQRLDN DVDLRGDQPS LGSFTYSTSA PGPALSPSVP LHYLPHDPLH QELSFGVPYS HMMPRRLSTQ RYRLQQPLPP PPPPPPPPY YPSFLPYFLS MLPMSPTAMG PTISLDLDVD DVEMENYEAL LNLAERLGDA KPRGLTKADI EQLPSYRFNP DSHQSEQTLC VVCFSDFEAR QLLRVLP CNH EFHTKCV DKW LKANRTCPIC RADASEVPRE AE
Source	Mammalian cell
Target Names	RNF44
Protein Names	Recommended name: RING finger protein 44
Expression Region	1-432
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This protein contains a RING finger, a motif present in a variety of functionally distinct proteins and known to be involved in protein-protein and protein-DNA interactions.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.