



Recombinant Human Ret finger protein-like 4B (RFPL4B)

Product Code	CSB-MP019603HU
Abbreviation	RFPL4B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6ZWI9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MAKRLQAELS CPVCLDFFSC SISLSCTHVF CFDCIQRYIL ENHDFRAMCP LCRDVVKVPA LEEWQVSVLT LMTKQHNSRL EQSLHVREEL RHFREDVTLD AATASSLLVF SNDLRSAQCK KIHHDLT KD P RLACVLGTPC FSSGQHYWEV EVGEVKSWSL GVCKEPADRK SNDLFPEHGF WISMKAGAIH ANTHLERIPA SPRLRRVGIF LDADLEEI QF FDVDNNVLIY THDGGFSLEL LCPFFCLELL GEGESGNVLT ICP
Source	Mammalian cell
Target Names	RFPL4B
Protein Names	Recommended name: Ret finger protein-like 4B Alternative name(s): RING finger protein 211
Expression Region	1-263
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.