



Recombinant Human Zinc finger protein 765 (ZNF765)

Product Code	CSB-BP027021HU
Abbreviation	ZNF765
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7L2R6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MALPQGLLTF RDVAIEFSQE EWKCLDPAQR TLYRDVMLEN YRNLVSLDIS SKCMMKEFSS TAQGNREVFH AGTSQRHESH HNGDFCFQDI DKDIHDIEFQ WQEDERNGHE ALMTKIKKLT GSTERYDQNY AGNKPVKYQL GFSFHSHP LHIFHTEEKI DNQVVKSIHD ASLVSTAQRI SCRPEETHISN DYGNNFLNSS LFTQKQEVHM REKSFQCNDG GKAYNCSSLL RKHQLIHLGE KQYKCDICGK VFNSKRYVAR HRRCHTGEKP YKCNECGKTF SQTYYLTCHR RLHTGEKPYK CEECDKAFHF KSKLQIHRRI HTGEKPYKCN ECGKTFSQKS YLTCHRRRLHT GEKPYKCNCGKTF SRKSHF TCHHRVHTGE KPYKCNCGKTF SHKSSLTY HRRRLHTEEKP YKCNECGKTF NQQLTLNCR LHSGEKPYK EEC DKAYSFK SNLEIHQKIH TEENPYKCN CGKTF SRSTSS LTYHHRRLHTG QKPYKCEDCD EAFSFKSNLE RHRRIYTGEK LHV
Source	Baculovirus
Target Names	ZNF765
Protein Names	Recommended name: Zinc finger protein 765
Expression Region	1-523
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.