



Recombinant Human Zinc finger protein 627 (ZNF627)

Product Code	CSB-MP026905HU
Abbreviation	ZNF627
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7L945
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MDSVAFEDVA VNFTLEEWAL LDPSQKNLYR DVMRETFRNL ASVGKQWEDQ NIEDPFKIPR RNISHIPERL CESKEGGQGE ETFSQIPDGI LNKKTPGVKP CESSVCGEVG MGPSSLNRHI RDHTGREPNE YQEYGKKSYS RNQCGRALSY HRSPVRETR HPGGKPYDCK ECGETFISLV SIRRHMLTHR GGVPYKCKVC GKAFDYPRLF RIHERSHTGE KPYECKQCGK AFSCSSYIRI HERTHTGDKP YECKQCGKAF SCSKYIRIHE RTHTGEKPYE CKQCGKAFRC ASSVRSHERT HTGEKLFECK ECGKALTCLA SVRRHMIKHT GNGPYKCKVC GKAFDFPSSF RIHERTHTGE KPYDCKQCGK AFSCSSSFRK HERIHTGEKP YKCTKCGKAF SRSSYFRIHE RTHTGEKPYE CKQCGKAFSR STYFRVHEKI HTGEKPYENP NPNASVVPVL S
Source	Mammalian cell
Target Names	ZNF627
Protein Names	Recommended name: Zinc finger protein 627
Expression Region	1-461
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.