



Recombinant Pig Phosphoglycerate kinase 1 (PGK1)

Product Code	CSB-EP017856PI-B
Abbreviation	PGK1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7SIB7
Product Type	Recombinant Protein
Immunogen Species	Sus scrofa (Pig)
Purity	>85% (SDS-PAGE)
Sequence	SLSNKLTL D KLDVKGKRVV MRVDFNVPMK NNQITNNQRI KAAIPSIKFC LDNGAKSVVL MSHLGRPDGI PMPDKYSLEP VAVELKSLPG KDVLFLKDCV GPEVEKACAD PAAGSVILLE NLRFHVEEEG KGKDASGSKV KADPAKIEAF RASLSKLGDV YVNDAFGTAH RAHSSMVGVN LPKKAGGFLM KKELNYFAKA LESPERPFLA ILGGAKVADK IQLINNMLDK VNEMIIGGGM AFTFLKVLNN MEIGTSLFDE EGSKIVKDL M SKAEKNGVKI TLPVDFVTAD KFDENAKIGQ ATVASGIPAG WMGLDCGPES SKKYSEAVAR AKQIVWNGPV GVFEWEAFAQ GTKALMDEVV KATSRGCITI IGGGDTATCC AKWNTEDKVS HVSTGGGASL ELLEKVLPG VDALS NV
Source	E.coli
Target Names	PGK1
Protein Names	Recommended name: Phosphoglycerate kinase 1 EC= 2.7.2.3
Expression Region	2-417
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This protein is a glycolytic enzyme that catalyzes the conversion of 1,3-diphosphoglycerate to 3-phosphoglycerate. The encoded protein may also act as a cofactor for polymerase alpha. This gene lies on the X-chromosome, while a related pseudogene also has been found on the X-chromosome and another on chromosome 19.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.