



# Recombinant Danio rerio Nucleolar protein 10 (nol10), partial

<b>Product Code</b>	CSB-BP015918DIL
<b>Abbreviation</b>	nol10
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q802W4
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Danio rerio (Zebrafish) (Brachydanio rerio)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MQVSSVNNVK IYNLSHGKSL PEWLSDRKKR VLQKKDVDIQ RRIELIQDFE MPTVCTSIRV SRDGQYILAA GTYKPRVRCY DTYQLSLKFE RCLDSDVVTF DILSDDYSKL VFLHIDRYVE FHSQHGHHYK TRIPKFGRDF SYHSPSCDLY FVGASSEVFR LNLEQGRFLN SLQTDAAEMN VCDINPVHQL FAAGTLEGRV DCWDPRVRTR VAALDCALSS ITDNTEVEGL PSVSALKFND SLGLAVGTST GQILVYDLRS SRPLLVKDHY YGLPIKSLHF HNSLDLVLSA DSKIIKMWNK DNGKVFSSIE PQANINDVCL YPASGMLFTA NEDPKMNTFY IPALGPAPRW CSFLDNLTEE LEENPESTIY DDYKFVTRKD LESLGLAHLI GSPLLRAYMH GFFMDIRLYH KVKTMVNPFA YEEYRKDKIR QKIEESRAQR VQLKKLPKVN KELALKLMEE DTELTKKKK KKANVAGNLL MDDRFKVMFE NPDYQVDERS EEFRLLNPII SNVAERRRKS LLEEEQEQAEE EQEPEGRGS SEDDSSDEDD KGWVQEVREQ RRLLRMEERE RSYIQRQERR QQDRNTRLQS SDTHTQQSHG AQKPRFYQIK SGEEFRSFSV VSRKQKTHKA SLEERLQRLE KSDTLNLNDT AVGSKQLTFT LKKTEKQRFQ QQAERDHHEE RRRIRRSAGH LHSNRGGGRG RGGGRGRGGG RGRGGGRGRR PF
<b>Source</b>	Baculovirus
<b>Target Names</b>	nol10
<b>Protein Names</b>	Recommended name: Nucleolar protein 10
<b>Expression Region</b>	1-722
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Partial
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.