



Recombinant Danio rerio Nucleolar protein 10 (nol10), partial

Product Code	CSB-MP015918DIL
Abbreviation	nol10
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q802W4
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MQVSSVNNVK IYNLSHGKSL PEWLSDRKKR VLQKKDVDIQ RRIELIQDFE MPTVCTSIRV SRDGQYILAA GTYKPRVRCY DTYQLSLKFE RCLDSDVVTF DILSDDYSKL VFLHIDRYVE FHSQHGHHYK TRIPKFGRDF SYHSPSCDLY FVGASSEVFR LNLEQGRFLN SLQTDAAEMN VCDINPVHQL FAAGTLEGRV DCWDPRVRTR VAALDCALSS ITDNTEVEGL PSVSALKFND SLGLAVGTST GQILVYDLRS SRPLLVKDHY YGLPIKSLHF HNSLDLVLSA DSKIIKMWNK DNGKVFSSIE PQANINDVCL YPASGMLFTA NEDPKMNTFY IPALGPAPRW CSFLDNLTEE LEENPESTIY DDYKFVTRKD LESLGLAHLI GSPLLRAYMH GFFMDIRLYH KVKTMVNPFA YEEYRKDKIR QKIEESRAQR VQLKKLPKVN KELALKLMEE DTELTKKKK KKANVAGNLL MDDRFKVMFE NPDYQVDERS EEFRLLNPII SNVAERRRKS LLEEEQEQAEE EQEPEGRGS SEDDSSDEDD KGWVQEVREQ RLLLRMEERE RSYIQRQERR QQDRNTRLQS SDTHTQQSHG AQKPRFYQIK SGEEFRSFSV VSRKQKTHKA SLEERLQRLE KSDTLNLNDT AVGSKQLTFT LKKTEKQRFQ QQAERDHEEE RRRIRRSAGH LHSNRGGGRG RGGGRGRGGG RGRGGGRGRR PF
Source	Mammalian cell
Target Names	nol10
Protein Names	Recommended name: Nucleolar protein 10
Expression Region	1-722
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.