



Recombinant Human Zinc finger protein 329 (ZNF329)

Product Code	CSB-MP026676HU
Abbreviation	ZNF329
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q86UD4
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MRLKMTTRNF PEREVPCDVE VERFTREVPC LSSLGDGWDC ENQEGHLRQS ALTLEKPGTQ EAICEYPGFG EHLIASSDLP PSQRVLATNG FHAPDSNVSG LDCDPALPSY PKSYADKRTG DSDACGKGFN HSMEVIHGRN PVREKPYKYP ESVKSFNHFT SLGHQKIMKR GKKS YEGKNF ENIFTLSSSL NENQRNLPGE KQYRCTECGK CFKRNSSLVL HHRHTHTGEKP YTCNECGKSF SKNYNLIVHQ RIHTGEKPYE CSKCGKAFSD GSALTQHQRI HTGEKPYECL ECGKTFNRNS SLILHQRTH HTGEKPYRCNEC GKPFTDISHL TVHLRIHTGE KPYECSKCGK AFRDGSYLTQ HERTHTGEKP FECAECGKSF NRNSHLIVHQ KIHSGEKPYE CKECGKTFIE SAYLIRHQRI HTGEKPYGCN QCQKLFRNIA GLIRHQRTH HTGEKPYECNQC GKAFRDSSCL TKHQRIHTKE TPYQCPECGK SFKQNSHLAV HQRLHSREGP SRCPQCGKMF QKSSSLVRHQ RAHLGEQPME T
Source	Mammalian cell
Target Names	ZNF329
Protein Names	Recommended name: Zinc finger protein 329
Expression Region	1-541
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.