



Recombinant Human Cerebellar degeneration-related protein 2-like (CDR2L)

Product Code	CSB-MP005104HU
Abbreviation	CDR2L
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q86X02
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MRRAAGMEDF SAEEEESWYD QQDLEQDLHL AEELGKTLLE RNKELEGLSLQ QMYSTNEEQV QEIEYLTKQL DTLRHVNEQH AKVYEQLDLT ARDLELTNHR LVLESKAAQQ KIHGLTETIE RLQAQVEELQ AQVEQLRGLQ QLRVLREKRE RRRTIHTFPC LKELCTSPRC KDAFRLHSSS LELGPRPLEQ ENERLQTLVG ALRSQVSQER QRKERAEREY TAVLQEYSEL ERQLCEMEAC RLRVQELEAE LLELQQMKQA KTYLLGPDDH LAEALLAPLT QAPEADDPQP GRGDDLGAQD GVSSPAASPG HVVRKSCSDT ALNAIVAKDP ASRHAGNLT L HANSVRKRG M SILREVDEQY HALLEKYEEL LSKCRQH GAG VRHAGVQTSR PISR DSSWRD LRG GEEGQGE VKAGEKSL SQ HVEAVDKRLE QSQPEYKALF KEIFSRIQKT KADINATKVK THSSK
Source	Mammalian cell
Target Names	CDR2L
Protein Names	Recommended name: Cerebellar degeneration-related protein 2-like Alternative name(s): Paraneoplastic 62 kDa antigen
Expression Region	1-465
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.