



Recombinant Human Zinc finger protein 227 (ZNF227)

Product Code	CSB-MP026595HU
Abbreviation	ZNF227
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q86WZ6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MPSQNYDLPQ KKQEKMTEKFKQ EAVTFKDVAV VFSREELRLL DLTQRKLYRD VMVENFKNLV AVGHLPFQPD MVSQLEAEEK LWMMETETQR SSKHQNKMET LQKFALKYLS NQELSCWQIW KQVASELTRC LQ GKSSQLLQ GDSIQVSENE NNIMNPKGDS SIYIENQEFP FWRTQHSCGN TYLSESQIQS RGKQIDVKNN LQIHEDFMKK SPFHHEIKTD TEPKPCGNE YGKIISDGSN QKLPLGEKPH PCGECGRGFS YSPRLPLHPN VHTGEKCFSSQ SSSLRTHQRI HPGEKLNKCH ESGDCFNKSS FHSYQSNHTG EKSRYCDSCG KGFSSSTGLI IHYRTHTEK PYKCEECGKC FSQSSNFQCH QRVHTEEKPY KCEECGKGF WSVNLRVHQR VHRGEKPYKC EECGKGFTQA AHFHIHQRVH TGEKPYKCDV CGKGFHNSP LICHRRVHTG EKPYKCEACG KGFTTRNTDLH IHFRVHTGEK PYKCKECCGKG FSQASNLQVH QNVHTGEKRF KCETCGKGF QSSKLQTHQR VHTGEKPYRC DVCCKDFSYS SNLKLHQVIH TGEKPYKCEE CGKGFWSRNS LHAHQRVHSG EKPYKCEQCD KSFSQAIDFR VHQRVHTGEK PYKCGVCGKG FSQSSGLQSH QRVHTGEKPY KCDVCGKGF YSSQFIYHQR GHTGEKPYKC EECGKGFGRS LNLRHHQRVH TGEKPHICEE CGKAFSLPSN LRVHLGVHTR EKLFKCEECG KGFSQSARLE AHQRVHTGEK PYKCDICDKD FRHRSRLTYH QKVHTGKKL
Source	Mammalian cell
Target Names	ZNF227
Protein Names	Recommended name: Zinc finger protein 227
Expression Region	1-799
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.