



Recombinant Human Zinc finger protein 554 (ZNF554)

Product Code	CSB-MP026834HU
Abbreviation	ZNF554
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q86TJ5
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MVTCAHLGRR ARLPAAQPSA CPGTCFSQEE RMAAGYLPRW SQELVTFEDV SMDFSQEEWE LLEPAQKNLY REVMLENYRN VVSLEALKNQ CTDVGIKEGP LSPAQTSQVT SLSSWTGYLL FQPVASSHLE QREALWIEEK GTPQASCSDW MTVLRNQDST YKKVALQEEP ASGINMIKLI REDGGWKQLE DSHEDPQGLL SQKASLHVVA VPQEKATAWH GFGENGNLSP ALVLSQGSSK GNHLCGSELD ITSLASDSVL NHHQLGYADR RPCESNECGN AIRQNSHFIQ HGGKMFVYLE NGQSLNHGMA LTIHNKINTA EKPFECHQCG KVFNRRHSLS EHQRIHTGEK PYECQECGRA FTHSSTLTRH LRTHTGEKPY GCGECGKAFN RISSLTQHQR IHTGEKPYKC EDCGKSFCQS SYLILHKRTH TGEKPYECSE CGKAFSDRSS LNQHERHTTG ENPYECKQCG RAFSQRSSLV RHERHTHTGEK PYRCQECGKA FSQSSSLVTH QKTHSSQKTY KIIDCGKAFY QNRHLIGY
Source	Mammalian cell
Target Names	ZNF554
Protein Names	Recommended name: Zinc finger protein 554
Expression Region	1-538
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.