



Recombinant Human Rab9 effector protein with kelch motifs (RABEPK)

Product Code	CSB-BP019230HU
Abbreviation	RABEPK
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7Z6M1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MKQLPVLEPG DKPRKATWYT LTPVGDSPCA RVGHSCSYLP PVGNAKRGKV FIVGGANPNR SFSDVHTMDL GKHQWDLDTG KGLLPRIYEHA SFIPSCTPDR IWFVGGANQS GNRNCLQVLN PETRTWTTPE VTSPPPSPRT FHTSSAAIGN QLYVFGGGER GAQPVQDTKL HVFDANTLTW SQPETLGNPP SPRHGHVMVA AGTKLFIHGG LAGDRFYDDL HCIDISDMKW QKLNPTGAAP AGCAAHSAVA MGKHVYIFGG MTPAGALDTM YQYHTEEQHW TLLKFDTLLP PGRLDHSMCI IPWPVTCASE KEDSNLTLN HEAEKEDSAD KVMSHSGDSH EESQTATLLC LVFGGMNTEG EIYDDCIVTV VD
Source	Baculovirus
Target Names	RABEPK
Protein Names	Recommended name: Rab9 effector protein with kelch motifs Alternative name(s): 40 kDa Rab9 effector protein p40
Expression Region	1-372
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.