



# Recombinant Mouse Hyaluronan and proteoglycan link protein 3 (Hapln3)

<b>Product Code</b>	CSB-YP010132MO
<b>Abbreviation</b>	Hapln3
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q80WM5
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	LPF YNGFYYSNDL SGRTLGNNGYG EGLFNGVKLV VETTEESLFS HQGASVTLPC HYHYEPALAS PRHVRVKWWK LSENGAPEQD VLVVIGQRHR SFGDYQGRVQ LRQDKQQEVS LELRDLRLED SGRYRCEVID GLEDESGLVE LELRGVVFPY QSREGRYQLN FHEAQQACQE QGAMVATFEQ LFRAWEEGLD WCNAGWLQDA SVQYPIVLPR QPCGGLGLAP GVRSYGQRHH RLHRYDVFCF AAALKGRVYY LENPKLTL EAREACQEDG AQISTVGQLF AAWKFRGLDR CDAGWLADGS ARYPIVHPRL NCGPPEPGVR TFGFDPHTR YGVYCYVQH
<b>Source</b>	Yeast
<b>Target Names</b>	Hapln3
<b>Protein Names</b>	Recommended name: Hyaluronan and proteoglycan link protein 3
<b>Expression Region</b>	18-359
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	This gene belongs to the hyaluronan and proteoglycan binding link protein gene family. This protein may function in hyaluronic acid binding and cell adhesion.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.