



Recombinant Dog Secreted frizzled-related protein 2 (SFRP2)

| | |
|--------------------------|---|
| Product Code | CSB-YP021139DO |
| Abbreviation | SFRP2 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q863H1 |
| Product Type | Recombinant Protein |
| Immunogen Species | Canis lupus familiaris (Dog) (Canis familiaris) |
| Purity | >85% (SDS-PAGE) |
| Sequence | LFFGQP DFSYKRSNCK PIPANLQLCH GIEYQNMRLP NLLGHETMKE VLEQAGAWIP LVMKQCHPDT KKFLCSLFAP VCLDDLDETI QPCHSLCVQV KDRCAPVMSA FGFPWPDMLE CDRFPQDNDL CIPLASSDHL LPATEEAPKV CEACKNKNE DNDIMETLCK NDFALKIKVK EITYINRDTK IILETKSKTI YKLNQVSRD LKKSVLWLKD SLQCTCEEMN DINAPYLVMG QKLG GELVIT SVKRWQKGQR EFKRISRSIR KLQC |
| Source | Yeast |
| Target Names | SFRP2 |
| Protein Names | Recommended name: Secreted frizzled-related protein 2 Short name= sFRP-2 |
| Expression Region | 25-294 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full Length of Mature Protein |
| Target Details | This gene encodes a member of the SFRP family that contains a cysteine-rich domain homologous to the putative Wnt-binding site of Frizzled proteins. SFRPs act as soluble modulators of Wnt signaling. Methylation of this gene is a potential marker for the presence of colorectal cancer. |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life |



of lyophilized form is 12 months at -20°C/-80°C.