



Recombinant Human Zinc finger protein 761 (ZNF761)

Product Code	CSB-MP027018HU
Abbreviation	ZNF761
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q86XN6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MAFSQGLLTF RDVAIEFSQE EWKCLDPAQR TLYRDVMLEN YRNLVSLDIS SKCTMKEFLS TAQGNREVFH AGTLQIHESH HNGDFCYQDV DKDIHDYEFQ WQEDERNGHE APMTKIKKLT GSTERYDQSH ARNKPIKDQL GSSFHSHLPE MHIFQTEEKI DNQVVKSIHD ASLVSTAQRI SCRPKTHISN NHGNNFWNSS LLTQKQEVHM REKSFQCNES GKAFNYSSLL RKHQIIHLAD KYKCDVCGKL FNQKRNLACH RRCHTGENPY KCNECGKTF SQTSSLTCHRR LHTGEKPYKC EECDKAFHFK SILERHRIIH TEEKPYKCNE CGKTFRQKSI LTRHRLHTG EKPYKCNECG KTFSHKSSLT CHHRLHTGEK PYKCNECGKT FSHKSSLTCH RRLHTGEKPY KCEECDKAYS FRSNFEIHRK IHTEDNAYKC NECGKTFSRT SSLTCHRRRH TGEQPYKCEE CDKAFRFKSN LERHRIHTG EKPYKCNECG KTFSRKSYLT CHHRLHTGEK AYKCNECSKT FSWKSSLTCH RRLHSGEKPY KCKEKGKTFN QQLTLKRHR LHSGENPYKC ESDKAYSFK SNLEIHQKIH TEQNPYKCNE CGKTF SRTSS LTCHRRRLHTG EKPYKCEECD KAFRVKSNLE GHRRIHTGEK PYKCNECGKT FSRKSYFVCH HRLHTGEKPY KCNECGKNFS QKSSLICHHR LHTGEKPYKC NECGKTF S QK SNLTCHRRRLH TGEKQV
Source	Mammalian cell
Target Names	ZNF761
Protein Names	Recommended name: Zinc finger protein 761
Expression Region	1-746
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.