



Recombinant Mouse Proline-rich acidic protein 1 (Prap1)

Product Code	CSB-BP018628MO
Abbreviation	Prap1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q80XD8
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	APAHQVPVKT KGKHFVPEQE TEKVWDTRAL EPLEKDNQLG PLLPEPKQKP AAAEKRPDA MTWVETEDIL SHLRSPQGP ELDLDSIDHP MSDDVQDEEV PQSRPILYRQ VLQGPEEDLD HLAHSMEDS
Source	Baculovirus
Target Names	Prap1
Protein Names	Recommended name: Proline-rich acidic protein 1 Alternative name(s): Pregnancy-specific uterine protein Uterine-specific proline-rich acidic protein
Expression Region	21-149
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.