



Recombinant 6-phosphogluconolactonase (pgl)

Product Code	CSB-MP017861SZB
Abbreviation	pgl
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q83LP7
Product Type	Recombinant Protein
Immunogen Species	Shigella flexneri
Purity	>85% (SDS-PAGE)
Sequence	MKQTVYIASP ESQQIHVWNL NHEGALTLTQ VVDVPGQVQP MIVSPDKRYL YVGRPEFRV LAYRIAPDDG ALTFAAESVL PGSPHISTD HQGQFVFGS YNAGNVSVTR LEDGLPVGVV DVVEGLDGCH SANISPDNRT LWVPALKQDR ICLFTVSDDG HLVAQDPAEV TVEGAGPRH MVFHPNEQYA YCVNELNSSV DVWELKDPHG NIECVQTLDM MPENFSDTRW AADIHITPDG RHLIYACDRTA SLITVFSVSE DGSVLSKEGF QPTETQPRGF NVDYSGKYLI AAGQKSHHIS VYEIVGEQGL LHEKGRYAVG QGPMWVVVNA H
Source	Mammalian cell
Target Names	pgl
Protein Names	Recommended name: 6-phosphogluconolactonase Short name= 6-P-gluconolactonase EC= 3.1.1.31
Expression Region	1-331
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.