



Recombinant Spermidine synthase (speE)

Product Code	CSB-YP022668SZB
Abbreviation	speE
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q83MF0
Product Type	Recombinant Protein
Immunogen Species	Shigella flexneri
Purity	>85% (SDS-PAGE)
Sequence	AEKKQWHET LHDQFGQYFA VDNVLYHEKT DHQDLIFEN AAFGRVMALD GVVQTTERDE FIYHEMMTHV PLLAHGHAKH VLIIGGGDGA MLREVTRHKN VESITMVEID AGVVSFCRQY LPNHNAGSYD DPRFQLVIDD GVNFNQTSQ TFDVIISDCT DPIGPGESLF TSAFYEGCKR CLNPGGIFVA QNGVCFLQQE EAIDSHRKLK HYFSDVGFYQ AAIPYYGGI MTFAWATDND ALRHLSTEII QARFLASGLK CRYYNPAIHT AAFALPQYLQ DALASQPS
Source	Yeast
Target Names	speE
Protein Names	Recommended name: Spermidine synthase EC= 2.5.1.16 Alternative name(s): Putrescine aminopropyltransferase Short name= PAPT SPDSY
Expression Region	2-288
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.