



Recombinant Human PHD finger protein 13 (PHF13)

Product Code	CSB-YP017902HU
Abbreviation	PHF13
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q86YI8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MDSDSCAAAF HPEEYSPSCK RRRTVEDFNK FCTFVLAYAG YIPYPKEELP LRSSPSPANS TAGTIDSDGW DAGFSDIASS VPLPVSDRCF SHLQPTLLQR AKPSNFLDR KKTDKLLKKK KRKRRDSDAP GKEGYRGGLL KLEAADPYVE TPTSPTLQDI PQAPSDPCSG WSDTTPSSGS CATVSPDQVK EIKTEGKRTI VRQGKQVVFR DEDSTGNDED IMVDSDDDSW DLVTCFCMKP FAGRPMIECN ECHTWIHLSC AKIRKSNVPE VFVCQKCRDS KFDIRRSNRS RTGSRKLFLLD
Source	Yeast
Target Names	PHF13
Protein Names	Recommended name: PHD finger protein 13 Alternative name(s): Survival time-associated PHD finger protein in ovarian cancer 1 Short name= SPOC1
Expression Region	1-300
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.