



Recombinant Human Zinc finger protein 410 (ZNF410)

Product Code	CSB-YP026730HU
Abbreviation	ZNF410
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q86VK4
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MLSDELESKP ELLVQFVQNT SIPLGQGLVE SEAKDITCLS LLPVTEASEC SRLMLPDDTT NHSNSSKEVP SSAVLRSLRV NVGPDGEETR AQTQKSPF LSTSESSLL QDLQPSDSTS FILLNLTRAG LGSSAEHLVF VQDEAEDSGN DFLSSESTDS SIPWFLRVQE LAHDSLIAAT RAQLAKNAKT SNGENVHLG SGDGQSKDSG PLPQVEKKLK CTVEGCDRTF VWPAHFYHL KTHRNDRSFI CPAEGCGKSF YVLQRLKVHM RTHNGEKPFM CHESGCGKQF TTAGNLKNHR RIHTGKPFM CEAQGCGRSF AEYSSLRKHL VVHSGEKPHQ CQVCGKTFSSQ SGRNVHMRK HHLQLGAAGS QEQEQTAEPL MGSSLLEEAS VPSKNLVSMN SQPSLGGESL NLPNTNSILG VDDEVLAEGS PRSLSSVPDV THHLVMTQSG RQSYEVSFLT AVNPQELLNQ GDLTERRT
Source	Yeast
Target Names	ZNF410
Protein Names	Recommended name: Zinc finger protein 410 Alternative name(s): Another partner for ARF 1 Zinc finger protein APA-1
Expression Region	1-478
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.