



Recombinant Phosphoglycolate phosphatase (c4155)

Product Code	CSB-EP017873EGX-B
Abbreviation	c4155
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8FCW0
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MNKFEDIRGV AFDLDGTLVD SAPGLAAVD MALYALELPV AGEERVITWI GNGADILMER ALTWARQERA TLRKTMGKPP VDDDIPAE EQ VRILRKLFD R YYGEVAEEGT FLFPHVADTL GALQAKGLPL GLVTNKPTPF VAPLLEALDI AKYFSVVIGG DDVQNKKPHP DPLLLVAERM GIAPQQMLFV GDSRNDIQAA KAAGCPSVGL TYGYNYGEAI DLSQPDVIYQ SLNDLLPALG LPHSENQESK ND
Source	E.coli
Target Names	c4155
Protein Names	Recommended name: Phosphoglycolate phosphatase Short name= PGP Short name= PGPase EC= 3.1.3.18
Expression Region	1-252
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.