



Recombinant Human Zinc finger protein 595 (ZNF595)

Product Code	CSB-EP026875HU-B
Abbreviation	ZNF595
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8IYB9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MELVTFRDVA IEFSPPEWKC LDPAQQNLYR DVMLENYRNL VSLGFVISNP DLVTCLEQIK EPCNLKIHET AAKPPAICSP FSQDLSPVQG IEDSFHKLIL KRYEKCGHEN LQLRKGCKRV NECKVQKGVN NGVYQCLSTT QSKIFQCNTC VKVFSKFSNS NKHKIRHTGE KPFKCTECGR SFYMSHLTQH TGIHAGEKPY KCEKCGKAFN RSTLSKHKR IHTGEKPYTC EECGKAFRRS TVLNEHKKIH TGEKPYKCEE CGKAFTRSTT LNEHKKIHTG EKPYKCKEKG KAFRWSTSMN EHKNIHTGEK PYKCKEKGKA FRQSRSLNEH KNIHTGEKPY TCEKCGKAFN QSSSLIHRH IHSEQKLYKC EECGKAFTWS SSLNKHKRIH TGEKPYACEE CGKAFYRSSH LAKHKRIHTG EKPYTCEEKG KAFNQSSTLI LHKRIHSGQK PYKCEEKGKA FTRSTTLNEH KKIHTGEKPY KCEEKGKAFI WSASLNEHKN IHTGEKPYKC KECGKAFNQS SGLIHRHSIH SEQKLYKCEE CGKAFTRSTA LNEHKKIHSG EKPYKCKEKG KAYNLSSTLT KHKRIHTGEK PFTCEEKGKA FNWSSSLTKH KIIHTGEKSY KCEEKGKAFN RPSTLTVHKR IHTGKEHS
Source	E.coli
Target Names	ZNF595
Protein Names	Recommended name: Zinc finger protein 595
Expression Region	1-648
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.