



# Recombinant Human Zinc finger protein 595 (ZNF595)

<b>Product Code</b>	CSB-MP026875HU
<b>Abbreviation</b>	ZNF595
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8IYB9
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MELVTFRDVA IEFSPPEWKC LDPAQQNLYR DVMLENYRNL VSLGFVISNP DLVTCLEQIK EPCNLKIHET AAKPPAICSP FSQDLSPVQG IEDSFHKLIL KRYEKCGHEN LQLRKGCKRV NECKVQKGVN NGVYQCLSTT QSKIFQCNTC VKVFSKFSNS NKHKIRHTGE KPFKCTECGR SFYMSHLTQH TGIHAGEKPY KCEKCGKAFN RSTLSKHKR IHTGEKPYTC EECGKAFRRS TVLNEHKKIH TGEKPYKCEE CGKAFTRSTT LNEHKKIHTG EKPYKCKEKG KAFRWSTSMN EHKNIHTGEK PYKCKEKGKA FRQSRSLNEH KNIHTGEKPY TCEKCGKAFN QSSSLIHRH IHSEQKLYKC EECGKAFTWS SSLNKHKRIH TGEKPYACEE CGKAFYRSSH LAKHKRIHTG EKPYTCEEKG KAFNQSSTLI LHKRIHSGQK PYKCEEKGKA FTRSTTLNEH KKIHTGEKPY KCEEKGKAFI WSASLNEHKN IHTGEKPYKC KECGKAFNQS SGLIHRHSIH SEQKLYKCEE CGKAFTRSTA LNEHKKIHSG EKPYKCKEKG KAYNLSSTLT KHKRIHTGEK PFTCEEKGKA FNWSSSLTKH KIIHTGEKSY KCEEKGKAFN RPSTLTVHKR IHTGKEHS
<b>Source</b>	Mammalian cell
<b>Target Names</b>	ZNF595
<b>Protein Names</b>	Recommended name: Zinc finger protein 595
<b>Expression Region</b>	1-648
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.