



Recombinant Human Exocyst complex component 8 (EXOC8)

Product Code	CSB-BP007887HU
Abbreviation	EXOC8
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8IYI6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAMAMSDSGA SRLRRQLESG GFEARLYVKQ LSQQSDGDRD LQEHRQRIQA LAEETAQNLK RNVYQNYRQF IETAREISYL ESEMYQLSHL LTEQKSSLES IPLTLLPAAA AAGAAAASGG EEGVGGAGGR DHLRGQAGFF STPGGASRDG SGPGEEGKQR TLTTLLEKVE GCRHLLTPG QYLVYNGDLV EYDADHMAQL QRVHGFLLMND CLLVATWLPQ RRGMYRYNAL YSLDGLAVVN VKDNPPMKDM FKLLMFPEsr IFQAENAKIK REWLEVLEDT KRALSEKRRR EQEEAAAPRG PPQVTSKATN PFEDDEEEEP AVPEVEEEKV DLSMEWIQEL PEDLDVCIQA RDFEGAVDLL DKLNHYLEDK PSPPPVKELR AKVEERVRQL TEVLVFELSP DRSLRGGPKA TRRAVSQIR LGQCTKACEL FLRNRAAAVH TAIRQLRIEG ATLLYIHKLC HVFFTSLLET AREFEIDFAG TDSGCYSAFV VWARSAMGMF VDAFSKQVFD SKESLSTAAE CVKVAKEHCQ QLGDIGLDTL FIIHALLVKD IQGALHSYKE IIIEATKHRN SEEMWRRMNL MTPEALGKLG EEMKSCGVSF FEQYTGDDCW VNLSYTVVAF TKQTMGFLEE ALKLYFPELH MVLLESLVEI ILVAVQHVDY SLRCEQDPEK KAFIRQNASF LYETVLPVVE KRFEEGVGKP AKQLQDLRNA SRLIRVNPES TTSVV
Source	Baculovirus
Target Names	EXOC8
Protein Names	Recommended name: Exocyst complex component 8 Alternative name(s): Exocyst complex 84 kDa subunit
Expression Region	1-725
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.