



# Recombinant Human Exocyst complex component 8 (EXOC8)

<b>Product Code</b>	CSB-EP007887HU-B
<b>Abbreviation</b>	EXOC8
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8IYI6
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MAMAMSDSGA SRLRRQLESG GFEARLYVKQ LSQQSDGDRD LQEHRQRIQA LAEETAQNLK RNVYQNYRQF IETAREISYL ESEMYQLSHL LTEQKSSLES IPLTLLPAAA AAGAAAASGG EEGVGGAGGR DHLRGQAGFF STPGGASRDG SGPGEEGKQR TLTTLLEKVE GCRHLETPG QYLVYNGDLV EYDADHMAQL QRVHGFLMND CLLVATWLPQ RRGMYRYNAL YSLDGLAVVN VKDNPPMKDM FKLLMFPEsr IFQAENAKIK REWLEVLEDT KRALSEKRRR EQEEAAPRG PPQVTSKATN PFEDDEEEEP AVPEVEEEKV DLSMEWIQEL PEDLDVCIAQ RDFEGAVDLL DKLNHYLEDK PSPPPVKELR AKVEERVRL TEVLVFELSP DRSLRGGPKA TRRAVSQIR LGQCTKACEL FLRNRAAAVH TAIRQLRIEG ATLLYIHKLC HVFFTSLLET AREFEIDFAG TDSGCYSAFV VWARSAMGMF VDAFSKQVFD SKESLSTAAE CVKVAKEHCQ QLGDIGLDTL FIIHALLVKD IQGALHSYKE IIIEATKHRN SEEMWRRMNL MTPEALGKLG EEMKSCGVSF FEQYTGDDCW VNLSYTVVAF TKQTMGFLLE ALKLYFPELH MVLLESLVEI ILVAVQHVDY SLRCEQDPEK KAFIRQNASF LYETVLPVVE KRFEEGVGKP AKQLQDLRNA SRLIRVNPES TTSVV
<b>Source</b>	E.coli
<b>Target Names</b>	EXOC8
<b>Protein Names</b>	Recommended name: Exocyst complex component 8 Alternative name(s): Exocyst complex 84 kDa subunit
<b>Expression Region</b>	1-725
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.