



Recombinant Argininosuccinate lyase (argH)

Product Code	CSB-MP002213EGX
Abbreviation	argH
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8FB96
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MALWGGRFTQ AADQRFKQFN DSLRFDYRLA EQDIVGSSVAW SKALVTVGVL TAEFAQLEE ALNVLLEDVR ARPQQILESD AEDIHSWVEG KLIDKVGQLG KKLHTGRSRN DQVATDLKLW CKDTVSELLT ANRQLQSALV ETAQNNQDAV MPGYTHLQRA QPVTFAHWCL AYVEMLARDE SRLQDALKRL DVSPGLGCGAL AGTAYEIDRE QLAGWLG FAS ATRNSLDSVS DRDHVLELLS AAAIGMVHLS RFAEDLIFFN TGEAGFVELS DRVTSGSSLM PQKKNPDALE LIRGKCGRVQ GALTGMMMTL KGLPLAYNKD MQEDKEGLFD ALDTWLDCLH MAALVLDGIQ VKRPRCQEAQ QQGYANATEL ADYLVAKGVP FREAHHIVGE AVVEAIRQ GK ALEDLPLSEL QKFSQVIGED VYPILSLQSC LDKRAAKGGV SPQQVAQAIA FAQARLG
Source	Mammalian cell
Target Names	argH
Protein Names	Recommended name: Argininosuccinate lyase Short name= ASAL EC= 4.3.2.1 Alternative name(s): Arginosuccinase
Expression Region	1-457
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.