



Recombinant Saimiri boliviensis boliviensis Cathepsin S (CTSS)

Product Code	CSB-BP006204SWC
Abbreviation	CTSS
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8HY82
Product Type	Recombinant Protein
Immunogen Species	Saimiri boliviensis boliviensis (Bolivian squirrel monkey)
Purity	>85% (SDS-PAGE)
Sequence	LPDSVD WREKGCVTEV KYQGSCGACW AFSAVGALEA QLKLKTGKLV SLSAQNLVDC SEKYGNKGCN GGFMTAEFQY IIDNKGIDSE ASYPYKATDQ KCQYDSKYRA ATCSKYTELP YGREDVLKEA VANKGPVCVG VDASHPSFFL YRSGVYYDPA CTQKVNHGVL VIGYGD LNGK EYWLVKNSWG SNFGEQGYIR MARNKGNHCG IASYPYPEI
Source	Baculovirus
Target Names	CTSS
Protein Names	Recommended name: Cathepsin S EC= 3.4.22.27
Expression Region	115-330
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.