



Recombinant Human Zinc finger protein 679 (ZNF679)

Product Code	CSB-MP026945HU
Abbreviation	ZNF679
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8IYX0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAKRPGSPGS REMGLLTFRD VVIEFSLEEW QCLDHAQQNL YRDVMLENYR NLVSLGIAVS KPDLITCLEQ NKEPWNIKRN EMVTKHPVMC SHFTQDLPPE LGIKDSLQKV IPRRYGKSGH DNLQVKTCKS MGECEVQKGG CNEVNQCLST TQNKIFQTHK CVKVFGKFSN SNRHKTRHTG KKHFKCKKYG KSFCMVSQLH QHIIHTREN SYQCEECGKP FNCSSSTLSKH KRIHTGEKPY RCEECGKAFT WSSTLTKHRR IHTGEKPYTC EECGQAFSRS STLANHKRIH TGEKPYTCEE CGKAFSLSSS LTYHKRIHTG EKPYTCEECG KAFNCSSSTLK KHKIIHTGEK PYKCKEKGKA FAFSSTLNTH KRIHTGEEPY KCEECDKAFK WSSSLANHKS MHTGEKPYKC E
Source	Mammalian cell
Target Names	ZNF679
Protein Names	Recommended name: Zinc finger protein 679
Expression Region	1-411
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.