



Recombinant Human RING finger protein 113B (RNF113B)

Product Code	CSB-MP019820HU
Abbreviation	RNF113B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8IZP6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAAPPSPGRT ADQADQVCTF LFKKPGRKGA AGLRKRACD PEHGESSSSG DEGDTVAQPP RVAPRPRGLH SWQKAAHGDR RGEEAAPESL DVVYRSTRSA KVPGPEDMGA TADFEQDTEK EHHTPTILKC SQRVQEALRG REHDHIYRGI HSYLRYLKP DTSMGNSSSG MARKGPIRAP GHLRATVRWD YQPDICKDYK ETGFCGFGDS CKFLHDRSDY KLGWEIEREL EEGRYCICED ENHEVGSEEE EIPFRCFICR QAFQNPVVTK CRHYFCESCA LEHFRATPRC YICDQPTGGI FNPAKELMAK LQKLQAAEGK KR
Source	Mammalian cell
Target Names	RNF113B
Protein Names	Recommended name: RING finger protein 113B Alternative name(s): Zinc finger protein 183-like 1
Expression Region	1-322
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.